



SIK2 blocking peptide (CDBP6108)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SIK2 salt-inducible kinase 2 [Homo sapiens (human)]
Official Symbol	SIK2
Synonyms	SIK2; salt-inducible kinase 2; QIK; SNF1LK2; LOH11CR1I; serine/threonine-protein kinase SIK2; SIK-2; SNF1-like kinase 2; qin-induced kinase; salt-inducible protein kinase 2; salt-inducible serine/threonine kinase 2; serine/threonine-protein kinase SNF1-like kinase 2
Entrez Gene ID	23235
mRNA Refseq	NM_015191
Protein Refseq	NP_056006
UniProt ID	Q9H0K1

Pathway	LKB1 signaling events
Function	ATP binding; magnesium ion binding; protein binding; protein serine/threonine kinase activity