



SIK3 blocking peptide (CDBP6109)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SIK3 SIK family kinase 3 [Homo sapiens (human)]
Official Symbol	SIK3
Synonyms	SIK3; SIK family kinase 3; L19; QSK; SIK-3; serine/threonine-protein kinase SIK3; salt-inducible kinase 3; serine/threonine-protein kinase QSK
Entrez Gene ID	23387
mRNA Refseq	NM_001281748
Protein Refseq	NP_001268677
UniProt ID	Q9Y2K2
Pathway	LKB1 signaling events

Function

ATP binding; magnesium ion binding; protein binding; protein serine/threonine kinase activity
