



SIPA1L1 blocking peptide (CDBP6115)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SIPA1L1 signal-induced proliferation-associated 1 like 1 [Homo sapiens (human)]
Official Symbol	SIPA1L1
Synonyms	SIPA1L1; signal-induced proliferation-associated 1 like 1; E6TP1; signal-induced proliferation-associated 1-like protein 1; SIPA1-like protein 1; high-risk human papilloma viruses E6 oncoproteins targeted protein 1
Entrez Gene ID	26037
mRNA Refseq	NM_001284245
Protein Refseq	NP_001271174
UniProt ID	O43166

Pathway	Rap1 signaling pathway
Function	GTPase activator activity; ephrin receptor binding; molecular_function