



SIPA1L2 blocking peptide (CDBP6116)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SIPA1L2 signal-induced proliferation-associated 1 like 2 [Homo sapiens (human)]
Official Symbol	SIPA1L2
Synonyms	SIPA1L2; signal-induced proliferation-associated 1 like 2; SPAL2; signal-induced proliferation-associated 1-like protein 2; SPA-1-like 2; SIPA1-like protein 2
Entrez Gene ID	57568
mRNA Refseq	NM_020808
Protein Refseq	NP_065859
UniProt ID	Q9P2F8
Pathway	Rap1 signaling pathway

Function	GTPase activator activity
----------	---------------------------
