



SIPA1L3 blocking peptide (CDBP6117)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SIPA1L3 signal-induced proliferation-associated 1 like 3 [Homo sapiens (human)]
Official Symbol	SIPA1L3
Synonyms	SIPA1L3; signal-induced proliferation-associated 1 like 3; SPAL3; signal-induced proliferation-associated 1-like protein 3; SPA-1-like 3; SIPA1-like protein 3; SPA-1-like protein 3
Entrez Gene ID	23094
mRNA Refseq	NM_015073
Protein Refseq	NP_055888
UniProt ID	O60292
Pathway	Rap1 signaling pathway

Function	GTPase activator activity
----------	---------------------------
