



SKA1 blocking peptide (CDBP6127)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SKA1 spindle and kinetochore associated complex subunit 1 [Homo sapiens (human)]
Official Symbol	SKA1
Synonyms	SKA1; spindle and kinetochore associated complex subunit 1; C18orf24; spindle and kinetochore-associated protein 1; spindle and KT (kinetochore) associated 1
Entrez Gene ID	220134
mRNA Refseq	NM_001039535
Protein Refseq	NP_001034624
UniProt ID	Q96BD8
Pathway	Cell Cycle; M Phase; Mitotic Anaphase; Mitotic Metaphase and Anaphase; Mitotic

Prometaphase; Resolution of Sister Chromatid Cohesion; Separation of Sister Chromatids

Function

microtubule binding; protein binding
