



SKA3 blocking peptide (CDBP6129)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a component of the spindle and kinetochore-associated protein complex that regulates microtubule attachment to the kinetochores during mitosis. The encoded protein localizes to the outer kinetochore and may be required for normal chromosome segregation and cell division. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SKA3 spindle and kinetochore associated complex subunit 3 [Homo sapiens (human)]
Official Symbol	SKA3
Synonyms	SKA3; spindle and kinetochore associated complex subunit 3; RAMA1; C13orf3; spindle and kinetochore-associated protein 3
Entrez Gene ID	221150
mRNA Refseq	NM_001166017

Protein Refseq	NP_001159489
UniProt ID	Q8IX90
Function	protein binding