



# SLAIN1 blocking peptide (CDBP6132)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|               |                                                            |
|---------------|------------------------------------------------------------|
| Conjugate     | Unconjugated                                               |
| Applications  | Used as a blocking peptide in immunoblotting applications. |
| Format        | Liquid                                                     |
| Concentration | 200 µg/mL                                                  |
| Size          | 0.05 mg                                                    |
| Preservative  | None                                                       |
| Storage       | -20°C                                                      |

## GENE INFORMATION

|                 |                                                                                  |
|-----------------|----------------------------------------------------------------------------------|
| Gene Name       | <a href="#">SLAIN1 SLAIN motif family, member 1 [ Homo sapiens (human) ]</a>     |
| Official Symbol | SLAIN1                                                                           |
| Synonyms        | SLAIN1; SLAIN motif family, member 1; C13orf32; SLAIN motif-containing protein 1 |
| Entrez Gene ID  | <a href="#">122060</a>                                                           |
| mRNA Refseq     | <a href="#">NM_001040153</a>                                                     |
| Protein Refseq  | <a href="#">NP_001035243</a>                                                     |
| UniProt ID      | Q8ND83                                                                           |