



SLAIN1 blocking peptide (CDBP6133)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SLAIN1 SLAIN motif family, member 1 [Homo sapiens (human)]
Official Symbol	SLAIN1
Synonyms	SLAIN1; SLAIN motif family, member 1; C13orf32; SLAIN motif-containing protein 1
Entrez Gene ID	122060
mRNA Refseq	NM_001040153
Protein Refseq	NP_001035243
UniProt ID	Q8ND83