



SLAIN2 blocking peptide (CDBP6134)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SLAIN2 SLAIN motif family, member 2 [Homo sapiens (human)]
Official Symbol	SLAIN2
Synonyms	SLAIN2; SLAIN motif family, member 2; KIAA1458; SLAIN motif-containing protein 2
Entrez Gene ID	57606
mRNA Refseq	NM_020846
Protein Refseq	NP_065897
UniProt ID	Q9P270
Function	protein binding