



SLAMF1 blocking peptide (CDBP6135)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SLAMF1 signaling lymphocytic activation molecule family member 1 [Homo sapiens (human)]
Official Symbol	SLAMF1
Synonyms	SLAMF1; signaling lymphocytic activation molecule family member 1; SLAM; CD150; CDw150; signaling lymphocytic activation molecule; IPO-3
Entrez Gene ID	6504
mRNA Refseq	NM_003037
Protein Refseq	NP_003028
UniProt ID	Q13291
Pathway	Measles

Function

antigen binding; protein binding; transmembrane signaling receptor activity
