



SLC22A17 blocking peptide (CDBP6137)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SLC22A17 solute carrier family 22, member 17 [Homo sapiens (human)]
Official Symbol	SLC22A17
Synonyms	SLC22A17; solute carrier family 22, member 17; BOCT; BOIT; 24p3R; NGALR; hBOIT; NGALR2; NGALR3; solute carrier family 22 member 17; 24p3 receptor; NGAL receptor; lipocalin-2 receptor; brain-type organic cation transporter; potent brain type organic ion transporter; neutrophil gelatinase-associated lipocalin receptor; solute carrier family 22 (organic cation transporter), member 17
Entrez Gene ID	51310
mRNA Refseq	NM_001289050
Protein Refseq	NP_001275979

UniProt ID	Q8WUG5
Function	transmembrane signaling receptor activity; transmembrane transporter activity