



SLC35D3 blocking peptide (CDBP6142)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SLC35D3 solute carrier family 35, member D3 [Homo sapiens (human)]
Official Symbol	SLC35D3
Synonyms	SLC35D3; solute carrier family 35, member D3; FRCL1; bA55K22.3; solute carrier family 35 member D3; frc, fringe-like 1; fringe connection-like protein 1
Entrez Gene ID	340146
mRNA Refseq	NM_001008783
Protein Refseq	NP_001008783
UniProt ID	Q5M8T2