



# SLC37A2 blocking peptide (CDBP6143)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

## GENE INFORMATION

Gene Name	<a href="#">SLC37A2 solute carrier family 37 (glucose-6-phosphate transporter), member 2 [ Homo sapiens (human) ]</a>
Official Symbol	SLC37A2
Synonyms	SLC37A2; solute carrier family 37 (glucose-6-phosphate transporter), member 2; pp11662; sugar phosphate exchanger 2; solute carrier family 37 member 2; solute carrier family 37 (glycerol-3-phosphate transporter), member 2
Entrez Gene ID	<a href="#">219855</a>
mRNA Refseq	<a href="#">NM_001145290</a>
Protein Refseq	<a href="#">NP_001138762</a>
UniProt ID	Q8TED4

Function	transporter activity
----------	----------------------

---