



SLFN12 blocking peptide (CDBP6144)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SLFN12 schlafen family member 12 [Homo sapiens (human)]
Official Symbol	SLFN12
Synonyms	SLFN12; schlafen family member 12; SLFN3
Entrez Gene ID	55106
mRNA Refseq	NM_001289009
Protein Refseq	NP_001275938
UniProt ID	Q8IYM2
Function	ATP binding