



SLFN14 blocking peptide (CDBP6145)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SLFN14 schlafen family member 14 [Homo sapiens (human)]
Official Symbol	SLFN14
Synonyms	SLFN14; schlafen family member 14
Entrez Gene ID	342618
mRNA Refseq	NM_001129820
Protein Refseq	NP_001123292
UniProt ID	P0C7P3
Function	ATP binding