



SLITRK2 blocking peptide (CDBP6149)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes an integral membrane protein that contains two N-terminal leucine-rich repeats domains and contains C-terminal regions similar to neurotrophin receptors. The encoded protein may play a role in modulating neurite activity. Alternatively spliced transcript variants encoding the same protein have been described.[provided by RefSeq, Feb 2010]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SLITRK2 SLIT and NTRK-like family, member 2 [Homo sapiens (human)]
Official Symbol	SLITRK2
Synonyms	SLITRK2; SLIT and NTRK-like family, member 2; CXorf2; SLITL1; SLIT and NTRK-like protein 2; slit-like 1; slit and trk like gene 2
Entrez Gene ID	84631
mRNA Refseq	NM_001144003

Protein Refseq [NP_001137475](#)

UniProt ID Q9H156
