



Slpi blocking peptide (CDBP6158)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Slpi secretory leukocyte peptidase inhibitor [Mus musculus (house mouse)]
Official Symbol	Slpi
Synonyms	Slpi; secretory leukocyte peptidase inhibitor; antileukoproteinase; ALP; secretory leukocyte protease inhibitor
Entrez Gene ID	20568
mRNA Refseq	NM_011414
Protein Refseq	NP_035544
UniProt ID	P97430
Function	enzyme binding; hydrolase activity; peptidase activity; peptidase inhibitor activity; serine-type

