



Diablo blocking peptide (CDBP6163)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Diablo diablo homolog (Drosophila) [Mus musculus (house mouse)]
Official Symbol	Diablo
Synonyms	Diablo; diablo homolog (Drosophila); Smac; AU040403; 0610041G12Rik; 1700006L01Rik; diablo homolog, mitochondrial; direct IAP binding protein with low PI; direct IAP-binding protein with low pl; second mitochondria-derived activator of caspase
Entrez Gene ID	66593
mRNA Refseq	NM_023232
Protein Refseq	NP_075721
UniProt ID	Q9JIQ3

Pathway

Apoptosis; Apoptotic factor-mediated response; Intrinsic Pathway for Apoptosis; Release of apoptotic factors from the mitochondria; SMAC binds to IAPs; SMAC-mediated apoptotic response; SMAC-mediated dissociation of IAP:caspase complexes; XPodNet - protein-protein interactions in the podocyte expanded by STRING
