



SOGA1 blocking peptide (CDBP6168)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SOGA1 suppressor of glucose, autophagy associated 1 [Homo sapiens (human)]
Official Symbol	SOGA1
Synonyms	SOGA1; suppressor of glucose, autophagy associated 1; SOGA; KIAA0889; C20orf117; protein SOGA1; SOGA family member 1; suppressor of glucose by autophagy; suppressor of glucose from autophagy; suppressor of glucose, autophagy-associated protein 1
Entrez Gene ID	140710
mRNA Refseq	NM_080627
Protein Refseq	NP_542194
UniProt ID	O94964

Functionmolecular_function
