



## SPATA3 blocking peptide (CDBP6173)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|               |                                                            |
|---------------|------------------------------------------------------------|
| Conjugate     | Unconjugated                                               |
| Applications  | Used as a blocking peptide in immunoblotting applications. |
| Format        | Liquid                                                     |
| Concentration | 200 µg/mL                                                  |
| Size          | 0.05 mg                                                    |
| Preservative  | None                                                       |
| Storage       | -20°C                                                      |

### GENE INFORMATION

|                 |                                                                                                                                                                                                                                                  |
|-----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene Name       | <a href="#">SPATA3 spermatogenesis associated 3 [ Homo sapiens (human) ]</a>                                                                                                                                                                     |
| Official Symbol | SPATA3                                                                                                                                                                                                                                           |
| Synonyms        | SPATA3; spermatogenesis associated 3; TSARG1; spermatogenesis-associated protein 3; testis spermatocyte apoptosis-related protein 1; testis and spermatogenesis cell-related protein 1; testis and spermatogenesis cell apoptosis related gene 1 |
| Entrez Gene ID  | <a href="#">130560</a>                                                                                                                                                                                                                           |
| mRNA Refseq     | <a href="#">NM_139073</a>                                                                                                                                                                                                                        |
| Protein Refseq  | <a href="#">NP_620712</a>                                                                                                                                                                                                                        |
| UniProt ID      | Q8NHX4                                                                                                                                                                                                                                           |

**Function**molecular\_function

---