



SPATA4 blocking peptide (CDBP6174)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SPATA4 spermatogenesis associated 4 [Homo sapiens (human)]
Official Symbol	SPATA4
Synonyms	SPATA4; spermatogenesis associated 4; SPEF1B; TSARG2; spermatogenesis-associated protein 4; testis and spermatogenesis cell related protein 2; testis and spermatogenesis cell-related protein 2; testis spermatocyte apoptosis-related gene 2 protein
Entrez Gene ID	132851
mRNA Refseq	NM_144644
Protein Refseq	NP_653245
UniProt ID	Q8NEY3

Function

molecular_function
