



SPATA6 blocking peptide (CDBP6175)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SPATA6 spermatogenesis associated 6 [Homo sapiens (human)]
Official Symbol	SPATA6
Synonyms	SPATA6; spermatogenesis associated 6; HASH; SRF1; SRF-1; spermatogenesis-associated protein 6; spermatogenesis-related factor-1
Entrez Gene ID	54558
mRNA Refseq	NM_001286238
Protein Refseq	NP_001273167
UniProt ID	Q9NWH7