



SPINK7 blocking peptide (CDBP6179)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SPINK7 serine peptidase inhibitor, Kazal type 7 (putative) [Homo sapiens (human)]
Official Symbol	SPINK7
Synonyms	SPINK7; serine peptidase inhibitor, Kazal type 7 (putative); ECG2; ECRG2; serine protease inhibitor Kazal-type 7; ECRG-2; esophagus cancer related gene 2; esophagus cancer-related gene-2; esophagus cancer-related gene 2 protein
Entrez Gene ID	84651
mRNA Refseq	NM_032566
Protein Refseq	NP_115955
UniProt ID	P58062

Function

protein binding; serine-type endopeptidase inhibitor activity
