



SPNS1 blocking peptide (CDBP6180)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SPNS1 spinster homolog 1 (Drosophila) [Homo sapiens (human)]
Official Symbol	SPNS1
Synonyms	SPNS1; spinster homolog 1 (Drosophila); LAT; nrs; SPIN1; SPINL; HSpin1; PP2030; protein spinster homolog 1; spinster-like protein 1
Entrez Gene ID	83985
mRNA Refseq	NM_001142448
Protein Refseq	NP_001135920
UniProt ID	Q9H2V7
Function	protein binding