



SPRYD4 blocking peptide (CDBP6186)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SPRYD4 SPRY domain containing 4 [Homo sapiens (human)]
Official Symbol	SPRYD4
Synonyms	SPRYD4; SPRY domain containing 4; SPRY domain-containing protein 4
Entrez Gene ID	283377
mRNA Refseq	NM_207344
Protein Refseq	NP_997227
UniProt ID	Q8WW59
Function	molecular_function