



TRIM51 blocking peptide (CDBP6187)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TRIM51 tripartite motif-containing 51 [Homo sapiens (human)]
Official Symbol	TRIM51
Synonyms	TRIM51; tripartite motif-containing 51; SPRYD5; TRIM51A; tripartite motif-containing protein 51; SPRY domain containing 5; SPRY domain-containing protein 5
Entrez Gene ID	84767
mRNA Refseq	NM_032681
Protein Refseq	NP_116070
UniProt ID	Q9BSJ1
Function	zinc ion binding