



SRPX blocking peptide (CDBP6190)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SRPX sushi-repeat containing protein, X-linked [Homo sapiens (human)]
Official Symbol	SRPX
Synonyms	SRPX; sushi-repeat containing protein, X-linked; DRS; ETX1; SRPX1; HEL-S-83p; sushi repeat-containing protein SRPX; sushi-repeat-containing protein, X-linked; sushi-repeat-containing protein, X chromosome; epididymis secretory sperm binding protein Li 83p
Entrez Gene ID	8406
mRNA Refseq	NM_001170750
Protein Refseq	NP_001164221
UniProt ID	P78539