



Il1rl1 blocking peptide (CDBP6192)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Il1rl1 interleukin 1 receptor-like 1 [Mus musculus (house mouse)]
Official Symbol	Il1rl1
Synonyms	Il1rl1; interleukin 1 receptor-like 1; T1; St2; DER4; Ly84; ST2L; Fit-1; T1/ST2; St2-rs1; interleukin-1 receptor-like 1; lymphocyte antigen 84
Entrez Gene ID	17082
mRNA Refseq	NM_001025602
Protein Refseq	NP_001020773
UniProt ID	P14719
Function	interleukin-1 receptor activity; interleukin-33 binding; interleukin-33 receptor activity