



STEAP3 blocking peptide (CDBP6199)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	STEAP3 STEAP family member 3, metalloredutase [Homo sapiens (human)]
Official Symbol	STEAP3
Synonyms	STEAP3; STEAP family member 3, metalloredutase; STMP3; TSAP6; AHMIO2; dudlin-2; metalloredutase STEAP3; pHyde; hTSAP6; hpHyde; dudlin 2; dudulin 2; dudulin-2; 1010001D01Rik; tumor suppressor pHyde; six transmembrane prostate protein 3; tumor suppressor-activated pathway protein 6; six-transmembrane epithelial antigen of prostate 3
Entrez Gene ID	55240
mRNA Refseq	NM_001008410
Protein Refseq	NP_001008410
UniProt ID	Q658P3

Pathway	Direct p53 effectors; Iron metabolism in placenta; Iron uptake and transport; Transferrin endocytosis and recycling; Transmembrane transport of small molecules; p53 signaling pathway
Function	ferric-chelate reductase activity; metal ion binding; protein binding