



## STOX1 blocking peptide (CDBP6207)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	The protein encoded by this gene may function as a DNA binding protein. Mutations in this gene are associated with pre-eclampsia/eclampsia 4 (PEE4). Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Used as a blocking peptide in immunoblotting applications.
<b>Format</b>	Liquid
<b>Concentration</b>	200 µg/mL
<b>Size</b>	0.05 mg
<b>Preservative</b>	None
<b>Storage</b>	-20°C

### GENE INFORMATION

<b>Gene Name</b>	<a href="#">STOX1 storkhead box 1 [Homo sapiens (human)]</a>
<b>Official Symbol</b>	STOX1
<b>Synonyms</b>	STOX1; storkhead box 1; PEE4; C10orf24; storkhead-box protein 1; winged-helix domain-containing protein
<b>Entrez Gene ID</b>	<a href="#">219736</a>
<b>mRNA Refseq</b>	<a href="#">NM_001130159</a>

<b>Protein Refseq</b>	<a href="#">NP_001123631</a>
<b>UniProt ID</b>	Q6ZVD7
<b>Function</b>	DNA binding