



STOX2 blocking peptide (CDBP6208)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	STOX2 storkhead box 2 [Homo sapiens (human)]
Official Symbol	STOX2
Synonyms	STOX2; storkhead box 2; storkhead-box protein 2
Entrez Gene ID	56977
mRNA Refseq	NM_020225
Protein Refseq	NP_064610
UniProt ID	Q9P2F5
Function	molecular_function