



SYNGR4 blocking peptide (CDBP6223)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes an integral membrane protein. The gene belongs to the synaptogyrin gene family. Like other members of the family the protein contains four transmembrane regions. The exact function of this protein is unclear. [provided by RefSeq, Jul 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SYNGR4 synaptogyrin 4 [Homo sapiens (human)]
Official Symbol	SYNGR4
Synonyms	SYNGR4; synaptogyrin 4; synaptogyrin-4
Entrez Gene ID	23546
mRNA Refseq	NM_012451
Protein Refseq	NP_036583

