



SYNPO blocking peptide (CDBP6224)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Synaptopodin is an actin-associated protein that may play a role in actin-based cell shape and motility. The name synaptopodin derives from the protein's associations with postsynaptic densities and dendritic spines and with renal podocytes (Mundel et al., 1997 [PubMed 9314539]).[supplied by OMIM, Mar 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SYNPO synaptopodin [Homo sapiens (human)]
Official Symbol	SYNPO
Synonyms	SYNPO; synaptopodin
Entrez Gene ID	11346
mRNA Refseq	NM_001109974
Protein Refseq	NP_001103444

UniProt ID	Q8N3V7
Function	actin binding; protein binding