



## SYNPO2 blocking peptide (CDBP6225)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

### GENE INFORMATION

Gene Name	<a href="#">SYNPO2 synaptopodin 2 [ Homo sapiens (human) ]</a>
Official Symbol	SYNPO2
Synonyms	SYNPO2; synaptopodin 2; synaptopodin-2; myopodin; genethonin-2
Entrez Gene ID	<a href="#">171024</a>
mRNA Refseq	<a href="#">NM_001128933</a>
Protein Refseq	<a href="#">NP_001122405</a>
UniProt ID	Q9UMS6
Function	14-3-3 protein binding; actin binding; muscle alpha-actinin binding; protein binding