



SYNPO2L blocking peptide (CDBP6226)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SYNPO2L synaptopodin 2-like [Homo sapiens (human)]
Official Symbol	SYNPO2L
Synonyms	SYNPO2L; synaptopodin 2-like; synaptopodin 2-like protein
Entrez Gene ID	79933
mRNA Refseq	NM_001114133
Protein Refseq	NP_001107605
UniProt ID	Q9H987
Function	actin binding