



SYPL1 blocking peptide (CDBP6229)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SYPL1 synaptophysin-like 1 [Homo sapiens (human)]
Official Symbol	SYPL1
Synonyms	SYPL1; synaptophysin-like 1; SYPL; H-SP1; synaptophysin-like protein 1; pantophysin
Entrez Gene ID	6856
mRNA Refseq	NM_006754
Protein Refseq	NP_006745
UniProt ID	Q16563
Function	transporter activity