



SYPL2 blocking peptide (CDBP6230)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SYPL2 synaptophysin-like 2 [Homo sapiens (human)]
Official Symbol	SYPL2
Synonyms	SYPL2; synaptophysin-like 2; MG29; synaptophysin-like protein 2; mitsugumin 29; mitsugumin-29
Entrez Gene ID	284612
mRNA Refseq	NM_001040709
Protein Refseq	NP_001035799
UniProt ID	Q5VXT5
Function	transporter activity