



TCTEX1D1 blocking peptide (CDBP6247)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TCTEX1D1 Tctex1 domain containing 1 [Homo sapiens (human)]
Official Symbol	TCTEX1D1
Synonyms	TCTEX1D1; Tctex1 domain containing 1; tctex1 domain-containing protein 1
Entrez Gene ID	200132
mRNA Refseq	NM_152665
Protein Refseq	NP_689878
UniProt ID	Q8N7M0