



TCTEX1D2 blocking peptide (CDBP6248)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TCTEX1D2 Tctex1 domain containing 2 [Homo sapiens (human)]
Official Symbol	TCTEX1D2
Synonyms	TCTEX1D2; Tctex1 domain containing 2; tctex1 domain-containing protein 2
Entrez Gene ID	255758
mRNA Refseq	NM_152773
Protein Refseq	NP_689986
UniProt ID	Q8WW35
Function	protein binding