



TCTE3 blocking peptide (CDBP6249)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TCTE3 t-complex-associated-testis-expressed 3 [Homo sapiens (human)]
Official Symbol	TCTE3
Synonyms	TCTE3; t-complex-associated-testis-expressed 3; TCTEX1D3; tctex1 domain-containing protein 3; tcte-3; Tctex1 domain containing 3; T-complex testis-specific protein 3; T-complex-associated testis-expressed protein 3
Entrez Gene ID	6991
mRNA Refseq	NM_174910
Protein Refseq	NP_777570
UniProt ID	Q8IZS6

Function	motor activity
----------	----------------
