



TCTEX1D4 blocking peptide (CDBP6250)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TCTEX1D4 Tctex1 domain containing 4 [Homo sapiens (human)]
Official Symbol	TCTEX1D4
Synonyms	TCTEX1D4; Tctex1 domain containing 4; tctex1 domain-containing protein 4; Tctex2 beta; tctex-2-beta; protein N22.1; novel Tctex-1 family domain-containing protein
Entrez Gene ID	343521
mRNA Refseq	NM_001013632
Protein Refseq	NP_001013654
UniProt ID	Q5JR98
Function	protein binding