



TCTN1 blocking peptide (CDBP6251)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a member of the tectonic family of secreted and transmembrane proteins. The orthologous gene in mouse is required for formation of most ventral cell types. It functions downstream of smoothened and rab23 to modulate hedgehog signal transduction. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TCTN1 tectonic family member 1 [Homo sapiens (human)]
Official Symbol	TCTN1
Synonyms	TCTN1; tectonic family member 1; TECT1; JBTS13; tectonic-1
Entrez Gene ID	79600
mRNA Refseq	NM_001082537

Protein Refseq [NP_001076006](#)

UniProt ID Q2MV58
