



TCTN2 blocking peptide (CDBP6252)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a type I membrane protein that belongs to the tectonic family. Studies in mice suggest that this protein may be involved in hedgehog signaling, and essential for ciliogenesis. Mutations in this gene are associated with Meckel syndrome type 8. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TCTN2 tectonic family member 2 [Homo sapiens (human)]
Official Symbol	TCTN2
Synonyms	TCTN2; tectonic family member 2; MKS8; TECT2; C12orf38; tectonic-2
Entrez Gene ID	79867
mRNA Refseq	NM_001143850

Protein Refseq [NP_001137322](#)

UniProt ID Q96GX1
