



TDRG1 blocking peptide (CDBP6255)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TDRG1 testis development related 1 (non-protein coding) [Homo sapiens (human)]
Official Symbol	TDRG1
Synonyms	TDRG1; testis development related 1 (non-protein coding); LINC00532
Entrez Gene ID	732253