



TEKT5 blocking peptide (CDBP6256)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TEKT5 tektin 5 [Homo sapiens (human)]
Official Symbol	TEKT5
Synonyms	TEKT5; tektin 5; CT149; tektin-5
Entrez Gene ID	146279
mRNA Refseq	NM_144674
Protein Refseq	NP_653275
UniProt ID	Q96M29