



CD248 blocking peptide (CDBP6258)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CD248 CD248 molecule, endosialin [Homo sapiens (human)]
Official Symbol	CD248
Synonyms	CD248; CD248 molecule, endosialin; TEM1; CD164L1; endosialin; 2610111G01Rik; CD164 sialomucin-like 1; CD248 antigen, endosialin; tumor endothelial marker 1
Entrez Gene ID	57124
mRNA Refseq	NM_020404
Protein Refseq	NP_065137
UniProt ID	Q9HCU0
Function	calcium ion binding; carbohydrate binding; extracellular matrix binding; molecular_function