



ARHGEF17 blocking peptide (CDBP6260)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ARHGEF17 Rho guanine nucleotide exchange factor (GEF) 17 [Homo sapiens (human)]
Official Symbol	ARHGEF17
Synonyms	ARHGEF17; Rho guanine nucleotide exchange factor (GEF) 17; TEM4; P164RHOGEF; p164-RhoGEF; rho guanine nucleotide exchange factor 17; tumor endothelial marker 4; 164 kDa Rho-specific guanine-nucleotide exchange factor; Rho-specific guanine-nucleotide exchange factor 164 kDa
Entrez Gene ID	9828
mRNA Refseq	NM_014786
Protein Refseq	NP_055601
UniProt ID	Q96PE2

Pathway	Cell death signalling via NRAGE; G alpha (12/13) signalling events; GPCR downstream signaling; NRAGE signals death through JNK; Regulation of RhoA activity; Rho GTPase cycle; Signal Transduction; Signaling by GPCR
Function	Rho guanyl-nucleotide exchange factor activity; guanyl-nucleotide exchange factor activity