



GPR124 blocking peptide (CDBP6261)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	GPR124 G protein-coupled receptor 124 [Homo sapiens (human)]
Official Symbol	GPR124
Synonyms	GPR124; G protein-coupled receptor 124; TEM5; G-protein coupled receptor 124; tumor endothelial marker 5
Entrez Gene ID	25960
mRNA Refseq	NM_032777
Protein Refseq	NP_116166
UniProt ID	Q96PE1
Pathway	Integrins in angiogenesis

FunctionG-protein coupled receptor activity
