



PLXDC1 blocking peptide (CDBP6263)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PLXDC1 plexin domain containing 1 [Homo sapiens (human)]
Official Symbol	PLXDC1
Synonyms	PLXDC1; plexin domain containing 1; TEM3; TEM7; plexin domain-containing protein 1; 2410003107Rik; tumor endothelial marker 3; tumor endothelial marker 7
Entrez Gene ID	57125
mRNA Refseq	NM_020405
Protein Refseq	NP_065138
UniProt ID	Q8IUK5
Function	protein binding; receptor activity